

60-100-0503



PRS-801 Resistance System Set

The Prostat PRS-801 is a unique, constant voltage, wide range ohmmeter with data logging, calculating and computer communication capabilities. Fully portable and battery operated, the PRS-801 measures resistance from <0.1 to 2.0×10^{14} ohms with measurement accuracy of $\pm 5\%$. It can be operated in Automatic, Automatic-Manual and Manual modes.

This sophisticated instrument digitally displays measurements in Ω , $K\Omega$, $M\Omega$, $G\Omega$ and $T\Omega$, and in exponential format, e.g., $1.5E08$.

It simultaneously displays data with a fully synchronous analog scale using $\times 1$, $\times 10$ and $\times 100$ multipliers, and a bank of color programmable LED's ($<10^3$ to $>10^{14}$ Ohm). Multiple test voltages include 0.01-10V, 10V and 100V ranges, which may be selected automatically by the instrument to comply with industry measurement practices or manually by the operator.

The PRS-801 records and stores up to 80 measurements, and calculates Minimum, Maximum and Average of all stored data.

Data may be downloaded to standard computer systems via the instrument's RS-232 Output Port, serial cable and supplied software for ESD Program documentation in Windows[®] Excel[®] spread sheets.

Fully automatic, the PRS-801 controls test voltage, resistance ranging and electrification period to make the most accurate and repeatable measurements available to ESD Auditors and laboratory technicians.

Specifications

Range:

Resistance from < 0.1 ($1.0E-1$) Ohms to 200 Tera ohms ($2.0E+14$ ohms). Maximum Resistivity with S11.11 concentric ring $2.0E+15$ ohms/square

TEST VOLTAGES

Automatic mode Default:

< 0.01 to 10 Volts Variable $1.0E-1$ to $1.0E+4$ Ohms

Constant Voltage:

10 Volts $\pm < 0.2$ Volts $1.0E+4$ to $1.0E+6$ Ohms

At 10 & 100 Volts 100 Volts \pm < 2.0 Volts 1.0E+6 to 2.0E+14 Ohms

Manual Mode:

< 0.01 to 10 Volts Variable 1.0E-1 to 1.0E+5 Ohms
 10 Volts \pm < 0.2 Volts 1.0E+3 to 1.0E+9 Ohms
 100 Volts \pm < 2.0 Volts 2.0E+5 to 2.0E+14 Ohms

Accuracy:

Overall: \pm <5% at ambient conditions (at 23°C and 30% Rh).

Specific Range Tolerances:

1.0E-1 to 1.0E+1 ohms: \pm 5% Corrected for Test Lead Resistance

1.0E+1 to 1.0E+12 ohms: \pm 2.0% with 10-foot Test Leads

1.0E+12 to 2.0E+14 ohms: \pm 40% or \pm 0.25 Decade with Grounded, Shielded Leads

Display:

Multi-function 2-5/8" x 1-5/8" Liquid Crystal Display with 0.5" digit height. Displays 3-1/2 digits in Ohms, or 1.0EXX in exponential format. Ohms Display indicators: Ω , K Ω , M Ω , G Ω and T Ω . Includes 19-segment analog scale (1-10 with 0.5 indication) with x1, x10, & x100 multipliers. Number of Data Points in Memory (0 – 80). Automatic Electrification Time (seconds), or Time required to Manually obtain steady state measurement. Displays data HOLD, BATTERY status, MIN, MAX, AVG, REC and Test Voltage (<10, 10, or 100 V)

LED Indicators:

14 Color Programmable LED's from $<10^3$ to $>10^{14}$ ohms. Colors (RED, GREEN, YELLOW, or Blank/OFF) may be operator programmed in SETUP mode.

Timer:

Times measurements In Seconds up to 99 seconds (Displayed on LCD)

Memory:

Register stores up to 80 data points (MEM # Displayed after RESET)

RS-232 Output:

Digital format: Exponential Power followed by Integer.

Response & Electrification:

Response from > 0.1 to $< 1.0E+06$ Ohms: < 2.0 seconds
 Average Measurement Period from 0.1 ohms to 1.0E12 Ohms **2.5** Seconds. Calculated Electrification Period per ESD S11.11, **7.5** seconds 0.1 ohms to 10E+12 Ohms.
 Programmed Electrification 1.0E+12 to $>2.0E14$ Ohms: **15.0+** seconds

Power:	Two 9 Volts alkaline batteries. Nominal battery life 25 hours.
Dimensions (W x L x D):	4.0" x 6.0" x 2.0" (10 cm x 15 cm x 5 cm)
Weight:	22 oz with batteries (624g)
Open Circuit Current (I):	< 4 mA @ 100 Volts

® Windows and Excel are the Registered Trade Marks of Microsoft Corporation

©2007 Prostat Corporation. All rights reserved.
Prostat is the registered trademarks of Prostat Corporation in the United States and other countries.

1072 Tower Lane, Bensenville, IL 60106 USA
Phone +1 (630) 238-8883 Fax +1 (630) 238-9717
Web access: <http://www.prostatcorp.com>

Specifications subject to change without notice.
Printed in U.S.A. REV2: 6/20/07